

# NEEDLE MICROTIP/SAFE WITH SAFETY MECHANISM



## MAIN SPECIFICATIONS

- Stainless steel needle AISI 304, lubricated with silicone oil
- Safety cap in colour code
- Safety mechanism
- Adaptable to both luer and luer lock syringes
- Safety mechanism: yes
- Type of needle: Safety hypodermic needles

Code	Description	Selling unit
-	NEEDLE MICROTIP/SAFE 16G-1,6x38 MM WHITE	CX. 50
--	NEEDLE MICROTIP/SAFE 18G-1,2x38 MM PINK	CX. 50
---	NEEDLE MICROTIP/SAFE 19G-1,1x38 MM PINK	CX. 50
----	NEEDLE MICROTIP/SAFE 19G-1,1x25 MM PINK	CX. 50
-----	NEEDLE MICROTIP/SAFE 20G-0,9x38 MM YELLOW	CX. 50
-----	NEEDLE MICROTIP/SAFE 20G-0,9x25 MM YELLOW	CX. 50

-----	NEEDLE MICROTIP/SAFE 21G-0,8x38 MM GREEN	CX. 50
-----	NEEDLE MICROTIP/SAFE 21G-0,8x25 MM GREEN	CX. 50
-----	NEEDLE MICROTIP/SAFE 22G-0,7x32 MM BLACK	CX. 50
-----	NEEDLE MICROTIP/SAFE 23G-0,6x25 MM BLUE	CX. 50
-----	NEEDLE MICROTIP/SAFE 23G-0,65x30 MM BLUE	CX. 50
-----	NEEDLE MICROTIP/SAFE 25G-0,5x16 MM ORANGE	CX. 50
-----	NEEDLE MICROTIP/SAFE 26G-0,45x12,7 MM BROWN	CX. 50
-----	NEEDLE MICROTIP/SAFE 27G-0,4x12,7 MM GREY	CX. 50
-----	NEEDLE MICROTIP/SAFE 29G-0,33x12,7 MM RED	CX. 50
-----	NEEDLE MICROTIP/SAFE 30G-0,30x8 MM YELLOW	CX. 50