


CATHETER VENICATH/SAFE WITHOUT WINGS AND INJECTION POINT



MAIN SPECIFICATIONS

- AISI 304 stainless steel needle with triple bevelling
- PP needle protection cap
- Transparent PP needle hub
- Catheter in radiopaque Polyurethane (available also in Teflon) with tapered tip
- PP Luer connector, suitable also for Luer Lock connectors, with PP Luer Lock white cap
- Safety mechanism: yes
- Type of needle: Standard I.V. Cannula

Code	Description	Selling unit
-	CATHETER VENICATH/SAFE WITHOUT WINGS AND INJECTION POINT 14 45 MM ORANGE	CX. 50
--	CATHETER VENICATH/SAFE WITHOUT WINGS AND INJECTION POINT 16G 45 MM GREY	CX. 50
---	CATHETER VENICATH/SAFE WITHOUT WINGS AND INJECTION POINT 17G 45 MM WHITE	CX. 50
----	CATHETER VENICATH/SAFE WITHOUT WINGS AND INJECTION	CX. 50

POINT 18G 45 MM GREEN

-----	CATHETER VENICATH/SAFE WITHOUT WINGS AND INJECTION POINT 20G 32 MM PINK	CX. 50
-----	CATHETER VENICATH/SAFE WITHOUT WINGS AND INJECTION POINT 22G 25 MM BLUE	CX. 50
-----	CATHETER VENICATH/SAFE WITHOUT WINGS AND INJECTION POINT 24G 19 MM YELLOW	CX. 50