

CATHETER HEMOFLON WITH WINGS AND INJECTION POINT



MAIN SPECIFICATIONS

- PUR Catheter
- Z-type valve (if present)
- Tapered tip
- Transparent flashback chamber
- Safety mechanism: no
- Type of needle: Standard I.V. Cannula

Code	Description	Selling unit
-	CATHETER HEMOFLON WITH WINGS AND INJECTION POINT 16G 45 MM GREY	CX. 50
--	CATHETER HEMOFLON WITH WINGS AND INJECTION POINT 17G 45 MM WHITE	CX. 50
---	CATHETER HEMOFLON WITH WINGS AND INJECTION POINT 18G 45 MM GREEN	CX. 50
----	CATHETER HEMOFLON WITH WINGS AND INJECTION POINT 18G 32 MM GREEN	CX. 50

-----	CATHETER HEMOFLON WITH WINGS AND INJECTION POINT 20G 32 MM PINK	CX. 50
-----	CATHETER HEMOFLON WITH WINGS AND INJECTION POINT 22G 25 MM BLUE	CX. 50
-----	CATHETER HEMOFLON WITH WINGS AND INJECTION POINT 24G 19 MM YELLOW	CX. 50