

BODEDEX - LIQUID DETERGENT FOR WASHING INSTRUMENTS



MAIN SPECIFICATIONS

- Instrument cleaner - patented system strong cleaning power
- High material compatibility
- For the active removal of blood, proteins, secretions and fat. The instrument detergent also gently dissolves the persistent residues of x-ray contrast media and biofilms
- Suitable for rigid and flexible optical instruments, intensive care and anesthesia materials, as well as surgical instruments, hospital units, primary health care and laboratory instruments.

Code	Description	Selling unit
16046021	BODEDEX FORTE 5L	UN
16046006	BODEDEX FORTE 2L	CX. 4