

# KORSOLEX READY TO USE - DISINFECTION OF INSTRUMENTS



## MAIN SPECIFICATIONS

- Ready-to-use aldehyde-based high-level instrument disinfectant
- Effective inactivation of all microorganisms
- Excellent material compatibility
- The ready-to-use Korsolex® is particularly suitable for the disinfection of all types of sensitive and heat-resistant instruments such as endoscopes, laparoscopes and anesthesia equipment.
- The solution does not need to be pre-diluted.

Code	Description	Selling unit
980 102	KORSOLEX READY TO USE 5L	UN