

# KORSOLEX PAA - DISINFECTION OF INSTRUMENTS



## MAIN SPECIFICATIONS

- Peracetic acid-based disinfectant for sensitive and heat-resistant instruments
- Ready to use, easy to apply solution
- Aldehyde-free
- Broad spectrum of effect, including virucidal efficacy
- Elevado nível de desinfeção
- Areas of application: flexible endoscopes, anesthesia material and surgical instruments; Thermoplastics and elastomers

Code	Description	Selling unit
980 408	KORSOLEX PAA WITH ACTIVATOR 5L	CX.2
980 409	KORSOLEX PAA 100 TEST STRIPS	CX. 6