

REINFORCED SURGICAL GOWN, STERILIZED



MAIN SPECIFICATIONS

- Reinforced, sterilized surgical gown with two hand towels
- Sterilized by Ethylene Oxide
- Material: SMMS Nonwoven Fabric
- Colors: Blue
- Packing: M - 40 sets (2 x 20), L - 40 sets (2 x 20), XL - 36 sets (2 x 18), XXL - 36 sets
- Available sizes: M, L, XL and XXL

Code	Description	Selling unit
12001 - S	SURGICAL GOWN SIZE S	UN
12001 - M	SURGICAL GOWN SIZE M	UN
12001 - L	SURGICAL GOWN SIZE L	UN
12001 - XL	SURGICAL GOWN SIZE XL	UN
12001 - XXL	SURGICAL GOWN SIZE XXL	UN