

HIDRAULIC STRETCHER



MAIN SPECIFICATIONS

- Steel tube structure, epoxy paint
- Phenolic compact bed with 100mm high density mattress
- Hydraulic lifting system via pedal
- Centralized brake and a directional wheel
- Adjustable height
- 2 sections
- Folding side guards
- Bullet holder for oxygen
- Soro and basket support for laundry
- Size: 2100 x 710 x (500 to 990) mm
- Weight: 98 kg

Code	Description	Selling unit
MA. 1790	HIDRAULIC STRETCHER	UN