

# EXAMINATION COUCH 3 SECTIONS



## MAIN SPECIFICATIONS

- Steel tube structure, epoxy paint
- Nappa-lined foam bed, gas spring adjustable headboard
- 4 swivel castors, central locking system and roll carrier
- Adjustable height by electric motor driven by a foot control
- Size: 1930 x 680 x (520 to 890) mm
- Weight: 65 kg

Code	Description	Selling unit
24019017	EXAMINATION COUCH 3 SECTIONS WITH BREATHING HOLE	UN
24009042	EXAMINATION COUCH 3 SECTIONS WITHOUT BREATHING HOLE	UN