

# ELECTRICAL EXAMINATION COUCH 3 SECTIONS



## MAIN SPECIFICATIONS

- Steel tube structure, epoxy paint
- Foam-lined foam bed, headboard and leg section, adjustable by gas spring
- With roller door, 4 swivel castors and central locking system
- Adjustable height by electric motor driven by a foot control
- Size: 1850 x 600 x (550 to 850) mm
- Weight: 65 kg

Code	Description	Selling unit
24019091	ELECTRICAL EXAMINATION COUCH 3 SECTIONS	UN