

# ELECTRICAL CHAIR "PODOLOGY"



## MAIN SPECIFICATIONS

- Steel tube structure, epoxy paint
- Foam-lined foam bed
- Equipped with two articulated armrests, adjustable in height and with 360 degree rotation
- With gas spring adjustable back section, with armrests
- With two independent leg sections adjustable by rack
- Height adjustable by electric motor driven by a foot control
- Size: 1850 x 600 x (540 to 850) mm
- Weight: 56 kg

Code	Description	Selling unit
CD. 1385	ELECTRICAL CHAIR "PODOLOGY"	UN