

# GYNECOLOGICAL TABLE MULTIFUNCTION



## MAIN SPECIFICATIONS

- Steel sheet structure, epoxy paint
- Foam-lined foam bed, 5 drawers, built-in step
- 3 sections being the headboard articulation by rack and extensible leg section, stainless steel skirting board
- Includes stainless steel waste tray, 2 adjustable thighs, and roller holder
- Size: 1700 x 590 x 850 mm.
- Weight: 87 kg

Code	Description	Selling unit
DV. 1672	GYNECOLOGICAL TABLE MULTIFUNCTIONS	UN