

# WOODEN EXAMINATION COUCH WITH SHELV



## MAIN SPECIFICATIONS

- Structure in beech wood
- Foam-lined foam bed, adjustable multi-position headboard with rack
- With roller door and bottom shelf
- Dimensions: 1850 x 600 x 750 mm
- Weight: 55.8 kg

| Code     | Description  | Selling unit |
|----------|--|--------------|
| DV. 1661 | WOODEN EXAMINATION COUCH WITH SHELV AND WITHOUT BREATHING HOLE | UN           |
| DV. 1663 | WOODEN EXAMINATION COUCH WITH SHELV AND WITH BREATHING HOLE    | UN           |