

# EXAMINATION COUCH



## MAIN SPECIFICATIONS

- Examination couch with steel tube structure, epoxy paint
- Foam-lined foam bed with adjustable headrest in various positions
- Removable feet with 4 plugs
- With roller bearings

Code	Description	Selling unit
24019037	EXAMINATION COUCH WITH BREATHING HOLE 1850X600X750 MM	UN
DV 1612	EXAMINATION COUCH WITH BREATHING HOLE 1850X600X600 MM	UN
DV 1637	EXAMINATION COUCH WITHOUT BREATHING HOLE 1850X600X750 MM	UN
24019092	EXAMINATION COUCH WITHOUT BREATHING HOLE 1850X600X600 MM	UN