

# CABINET FOR MEDICINES WITH 1 DOOR



## MAIN SPECIFICATIONS

- Structure made of sheet steel with rounded edges, epoxy paint
- 1 door with key, side and 4 glass shelves adjustable in height
- Size: 360 x 520 x 1450 mm
- Weight: 35.6 kg

Code	Description	Selling unit
24019018	CABINET FOR MEDICINES WITH 1 DOOR	UN
AR 1101	CABINET FOR MEDICINES WITH 1 DOOR AND FOOTS IN STAINLESS STEEL	UN