

## DRESSING TABLE



### MAIN SPECIFICATIONS

- Structure in steel tube, epoxy or stainless steel, with shelf in steel plate
- Includes 2 stainless steel removable tubs
- With 4 wheels
- Dimensions: 420 x 620 x 850 mm

Code	Description	Selling unit
ME 1840	DRESSING TABLE EPOXY	UN
ME 1841	DRESSING TABLE STAINLESS STEEL	UN
ME 1845	DRESSING TABLE EPOXY AND GLASS TOP	UN