

# ARTICULATED MANUAL BED WITH SIDE RAILS, HEAD AND FOOT BOARDS IN BEECH WOOD



## MAIN SPECIFICATIONS

- Steel tube structure, epoxy paint
- Deck of 4 planes with 3 joints, bed in wooden rulers
- 4 wheels being 2 with brake
- With headboards and guards in beech wood
- Dimensions: 1960 x 1020 x 500 mm
- Weight: 75 kg

Code	Description	Selling unit
CM 6002	ARTICULATED MANUAL BED WITH SIDE RAILS, HEAD AND FOOT BOARDS IN BEECH WOOD	UN