

# ARTICULATED MANUAL BED IN BEECH WOOD WITH SIDE RAILS



## MAIN SPECIFICATIONS

- Outer structure in beech wood without wheels, square feet, steel tube platform, epoxy paint
- 4 planes with 3 joints, bed in wooden slats, headboards with opening and side guards in beech wood
- Dimensions: 2020 x 1020 x 500 mm
- Weight: 66 kg

Code	Description	Selling unit
CM 6066	ARTICULATED MANUAL BED IN BEECH WOOD WITH SIDE RAILS	UN