

DRUG-SCREEN - DRUG DETECTION TEST / SALIVA METHOD



MAIN SPECIFICATIONS

- Contents: Packing with 15 Drug-Screen Multi Saliva; Instructions for use
- Drug-Screen-Multi 5CC Saliva Drug Rapid Test is a highly accurate and reliable cassette test.
- Test time: Accurate results after 10 minutes
- Storage: 2 ° to 30 ° C
- Validity: 12-18 months
- Parameters (cut-off): AMP50 Amphetamines (50ng / ml); COC20 Cocaine (20ng / ml); MET50 Methamphetamines (50ng / ml); MOR40 Opiates / Morphine (40ng / ml); THC12 Cannabinoids (12ng / ml)
- Advantages:
 - Easy to handle and simple read results
 - Accurate and accurate results after 10 minutes
 - Detailed Product Description
 - 99% agreement with GC / MS reference method
 - Limits of detection according to the international standard of SAMHSA / NIDA

Code	Description	Selling unit
-	DRUG-SCREEN MULTI SALIVA 5 CC	CX.15