

ABSORBENT UNDERPADS NEXDERMA



MAIN SPECIFICATIONS

- Material: Cellulose pulp, polypropylene, polyethylene.
- Color: White absorbent area and blue waterproof polyethylene.
- Dimension: 60 x 60 cm, 60 x 90 cm and 80 x 180 cm.
- Sterility: Non-sterile.
- Features: Waterproof, softness and comfort with uniform absorption.
- Model: Side flaps that help fix and prevent leakage. Absorbent area soft and developed for comfort.
- Absorbent surface of high absorption and allows good distribution of fluids.
- Developed for use in urinary and / or fecal incontinence, keeping the skin of the server protected from moisture, avoiding skin lesions resulting from contact with moisture and providing comfort to the user.
- Keep in a dry and dust free place. Avoid moisture, direct sunlight and fluorescent light. Temperature of not less than 10 ° C and not more than 40 ° C.
- Shelf life: 5 years, after manufacturing date.

Code	Description	Selling unit
R6090	ABSORBENT UNDERPAD 60X90 CM	CX. 150
R80180	ABSORBENT UNDERPAD 80X180 CM	CX. 100

