

WINGED NEEDLE



MAIN SPECIFICATIONS

- Winged needle (also called butterflies) for short term venous access.
- Flexible wings for an easier fixing.
- Smooth insertion needle.
- Siliconized needle with micro-bevel for flexible wings for atraumatic insertion for matching the puncture site
- Gauges from 27G to 19G - Wings color-coded.
- Extender allows connection to a pump infusion systems.
- 30 cm long tube, thus distancing the insertion point manipulations.
- Extraflexível tube ensuring total immobility of the needle into the vein.
- The tube does not get wrinkles even with the use of tweezers or clamps.
- Individually packaged in sterile blister for easy opening.

Code	Description	Selling unit
07001035	AGULHA EPICRANEANA VENOFIX G-21 0,8X20 mm	Cx.50
07001036	AGULHA EPICRANEANA VENOFIX G-23 0,65x20 mm	Cx.50

07001031	AGULHA EPICRANEANA VENOFIX G-19 1,1x20 mm	CX.50
----------	---	-------

07001032	AGULHA EPICRANEANA VENOFIX G-25 0,5x15 mm	CX.50
----------	---	-------
