


# OXYGEN MASK MEDIUM CONCENTRATION



## MAIN SPECIFICATIONS

- Mask for oxygen administration with metal clip to better fit nose.
- Transparent cover for better visual display of the patient face and vital signs, the mask is made of colorless and transparent PVC to permit visual inspection.
- With elastic band for fixing the mask to the head.
- With 2 m tube.
- Single use.
- Swivel Connector - swivel connector fits all oxygen tubes that meet ISO standards and offers maximum comfort, to be able to turn around the required position.

Code	Description	Selling unit
39279001	ADULT OXYGEN MASK MEDIUM CONCENTRATION WITH TUBING	CX.100
39011047	PEDIATRIC OXYGEN MASK MEDIUM CONCENTRATION WITH TUBING	CX.100
39011021	ADULT OXYGEN MASK MEDIUM CONCENTRATION WITH 7 VENTURI ADAPTERS	CX.100

FACTORMED  
Produtos Médicos & Hospitalares