

# GAZE COMPRESSE NON STERILE



## MAIN SPECIFICATIONS

- Composition: 100% cotton wool, according to European Standard EN14079.
- Available in various dimensions, folds, layers, wire numbers / cm2.
- Wired or wireless detectable to the RX.
- Presentation in plastic bags.
- Packaging in plastic bag.

Code	Description	Selling unit
01043007	GAZE COMPRESSE NON STERILE 8 LAYERS 17 WIRE/CM2 5x5 CM	CX. 100
01043006	GAZE COMPRESSE NON STERILE 8 LAYERS 17 WIRE/CM2 7,5X7,5 CM	CX. 200
01043002	GAZE COMPRESSE NON STERILE 8 LAYERS 17 WIRE/CM2 10x10 CM	CX. 200
01043003	GAZE COMPRESSE NON STERILE 8 LAYERS 17 WIRE/CM2 10x20 CM	CX. 200
01043008	GAZE COMPRESSE NON STERILE 8 LAYERS 17 WIRES/CM2 20X15 CM	CX.100
01043005	GAZE COMPRESSE NON STERILE 8 LAYERS 17 WIRES/CM2 20X20 CM	CX. 200

01043012	GAZE COMPRESSE NON STERILE 8 LAYERS 17 WIRE/CM2 40X20 CM	CX.100
01043009	GAZE COMPRESSE NON STERILE 16 LAYERS 17 WIRE/CM2 10x10 CM	CX.100
01043051	GAZE COMPRESSE NON STERILE WITH RX CONTRAST 8 LAYERS 17 WIRES/CM2 40X20CM	CX.100
01043052	GAZE COMPRESSE NON STERILE WITH RX CONTRAST 8 LAYERS 17 WIRES/CM2 10X10CM	CX. 200